

REGULAR

Embossed paper examination rolls

- **1 ply** or **2ply**: 100% pure paper pulp
- **Glue laminated**
- Embossed
- Soft
- Maximum absorbency
- Perforated
- Wrapped separately with PE shrink film



Suggested use

MEDICAL

BEAUTY

HEALTH

TATTOO



white

	HEIGHT	LENGTH	COLOURS	PLIES	PERFORATED PER	PIECES	ROLLS/CARTON	CARTONS/PALLET
STANDARD	40 cm	50 m	● PRINT	2	38 cm	132	12	40
		80 m	● PRINT	2	38 cm	210	8	40
		100 m	● PRINT	2	38 cm	264	6	40
	50 cm	50 m	● PRINT	2	38 cm	132	12 or 6	32 or 64
		80 m	● PRINT	2	38 cm	210	8	32
		100 m	● PRINT	2	38 cm	264	6	32
	58 cm	50 m	● PRINT	2	38 cm	132	12 or 6	32 or 64
		80 m	● PRINT	2	38 cm	210	8	32
		100 m	● PRINT	2	38 cm	264	6	32
68 cm	50 m	● PRINT	2	38 cm	132	12	24	
	80 m	● PRINT	2	38 cm	210	8	24	
	100 m	● PRINT	2	38 cm	264	6	24	
NORMAL	40 cm	50 m	● PRINT	1	38 cm	132	12	40
		80 m	● PRINT	1	38 cm	210	8	40
		100 m	● PRINT	1	38 cm	264	6	40
	50 cm	50 m	● PRINT	1	38 cm	132	12 or 6	32 or 64
		80 m	● PRINT	1	38 cm	210	8	32
		100 m	● PRINT	1	38 cm	264	6	32
	58 cm	50 m	● PRINT	1	38 cm	132	12 or 6	32 or 64
		80 m	● PRINT	1	38 cm	210	8	32
		100 m	● PRINT	1	38 cm	264	6	32
	68 cm	50 m	● PRINT	1	38 cm	132	12	24
		80 m	● PRINT	1	38 cm	210	8	24
		100 m	● PRINT	1	38 cm	264	6	24

Custom width & length upon request.

Ability to print other colors or designs (e.g. logo) upon request.

Custom perforation upon request.

SERIES	TYPES	COMPOSITION	ICONS	PLIES	LAMINATION	SOFTNESS	ABSORBENCY	RESILIENCE	PE	WATERPROOF / HYDROPHOBIC	COLOURS	SUGGESTED USE
ELITE STANDARD	Examination roll	100% pure paper pulp & LDPE		2	Thermal	●●●○○	●●●●○	●●●●●	●	Waterproof		MEDICAL, HEALTH, BEAUTY, TATTOO
PREMIUM PLUS	Examination roll	100% pure paper pulp & HDPE		2	Glue	●●●○○	●●●●○	●●●●●	●	Waterproof		MEDICAL, HEALTH, BEAUTY, TATTOO
PREMIUM STANDARD	Examination roll	100% pure paper pulp & HDPE		2	Glue	●●●○○	●●●●○	●●●○○	●	Waterproof		MEDICAL, HEALTH, BEAUTY, TATTOO
ECO	Examination roll	Recycled Pulp & HDPE		2	Glue	●●○○○	●●●●○	●●○○○	●	Waterproof		MEDICAL, HEALTH, BEAUTY, TATTOO
REGULAR NORMAL	Examination roll	100% pure paper pulp		1	-	●●●○○	●●●●○	●●●○○	-	-		MEDICAL, HEALTH, BEAUTY, TATTOO
REGULAR STANDARD	Examination roll	100% pure paper pulp		2	Glue	●●●○○	●●●●○	●●○○○	-	-		MEDICAL, HEALTH, BEAUTY, TATTOO
CLASSIC NORMAL	Examination roll	100% pure paper pulp		1	-	●●●○○	●●●●○	●●●○○	-	-		MEDICAL, HEALTH, BEAUTY, TATTOO
CLASSIC STANDARD	Examination roll	100% pure paper pulp		2	-	●●●○○	●●●●○	●●○○○	-	-		MEDICAL, HEALTH, BEAUTY, TATTOO
SUPREME EXTRA PLUS	Examination roll	100% polypropylene		1	-	●●●●●	-	●●●●●	-	Hydrophobic		MEDICAL, HEALTH, BEAUTY, TATTOO
SUPREME PLUS	Examination roll	100% polypropylene		1	-	●●●●●	-	●●●●○	-	Hydrophobic		MEDICAL, HEALTH, BEAUTY, TATTOO
SUPREME STANDARD	Examination roll	100% polypropylene		1	-	●●●●●	-	●●●○○	-	Hydrophobic		MEDICAL, HEALTH, BEAUTY, TATTOO
PRESTIGE	Examination roll	100% polypropylene & LDPE		2	Thermal	●●●●●	-	●●●●●	●	Waterproof		MEDICAL, HEALTH, BEAUTY, TATTOO
HYPER	Examination roll	100% Airlaid & HDPE		2	-	●●●●○	●●●●●	●●●●●	●	Waterproof		MEDICAL, HEALTH, BEAUTY, TATTOO
HYPER	Examination roll	100% Airlaid		1	-	●●●●○	●●●●●	●●●●●	-	-		MEDICAL, HEALTH, BEAUTY, TATTOO
HYPER	Cleaning Wipes	100% Airlaid		1	-	●●●●○	●●●●●	●●●●●	-	-		MEDICAL, HEALTH, BEAUTY, TATTOO, HO.RE.CA., CLEANING, HOME USE
GENTLE Alfa-Alfa	Cellulose wadding paper	100% pure paper pulp		2	-	●●●●●	●●●○○	-	-	-		MEDICAL, HEALTH, CLEANING
GENTLE Extra-Soft	Cellulose wadding paper	recycled pulp		2	-	●●●●○	●●●○○	-	-	-		MEDICAL, HEALTH, CLEANING
GENTLE Alfa-Soft	Cellulose wadding paper	recycled pulp		2	-	●●●○○	●●●○○	-	-	-		MEDICAL, HEALTH, CLEANING
GENTLE Extra-Extra	Cellulose wadding paper	recycled pulp		2	-	●●●●○	●●●○○	-	-	-		MEDICAL, HEALTH, CLEANING
DELICATE	Industrial paper roll	100% pure paper pulp		2	-	●●●○○	●●●●○	●●●●○	-	-		HEALTH, HO.RE.CA., HOME USE, BEAUTY, TATTOO, CLEANING
DELICATE	Industrial paper wipes	100% pure paper pulp		2	-	●●●○○	●●●○○	●●●○○	-	-		HEALTH, HO.RE.CA., HOME USE, BEAUTY, TATTOO, CLEANING