



ResMed

Actual Size: 140 mm (Width)
x 175 mm (Height)

Ruler is for scaling purposes.

How to use the fitting template:

1. Please ensure the template is printed to scale
2. Cut out along dotted line
3. Refer to diagram on template to determine appropriate cushion size



AirTouch™ N20

AirFit™ N20

Nasal mask
Nasal mask **for Her**

© 2020 ResMed Pty Ltd. All Rights Reserved.
Specifications may change without notice.
Product availability may vary across regions.
10111885/1 2020-02

Rx Only 0123



Sample size check